

# Electrochemical analysis chip 7 + 8: chip level measurement

The measuring adapter can directly accommodate and contact an electrochemical analytical chip

Two new layouts have been implemented at Fraunhofer IPMS for electrochemical analysis. The chip size is  $5 \times 5 \text{ mm}^2$  in each case. The chips have well-defined gold electrodes and one silver electrode each, which can be subsequently chlorinated to be used as a reference electrode. The chips are designed to be used with the Fraunhofer IPMS measurement adapter (see above), which contacts the electrodes via gold pins and can be connected directly to a potentiostat.

## Chip description electrochemical analytical chip 7

The analytical chip 7 consists of 4 gold electrodes and a silver electrode arranged in the center (see above). The electrodes have a diameter of approx. 107  $\mu$ m. Each electrode can be connected individually via a 700 x 700  $\mu$ m contact pad.



Chip layout and bond pad grid, valid for chips 7 and 8



Microscope image of the electrochemical analytical chip 7

#### Contact

Dr. Olaf Rüdiger Hild +49 351 8823-450 olaf.hild@ipms.fraunhofer.de

Fraunhofer Institute for Photonic Microsystems IPMS Maria-Reiche-Strasse 2 01109 Dresden Germany

www.ipms.fraunhofer.de



SEM images of a micro gold electrode (left) and a micro silver electrode (center) Right: SEM image of the electrode structures with the silver electrode in the center surrounded by 4 gold electrodes

#### Chip description electrochemical analytical chip 8

Analytical chip 8, shown above left and center, has a threeelectrode arrangement with a centered small silver electrode (1 960  $\mu$ m<sup>2</sup>) that can be used as a reference electrode, a gold working electrode (181 000  $\mu$ m<sup>2</sup>), and a counter electrode (181 000  $\mu$ m<sup>2</sup>) also made of gold.



Microscope image of the electrochemical analytical chip 8

### **Contact pads of the chips**

With a size of  $0.7 \times 0.7 \text{ mm}^2$  each, the contact pads of chip 7 and 8 are large enough to be contacted with a needle adapter or even small, back-insulated crocodile clips.



SEM image of a gold contact pad, on the left in the overview, right a section of a pad corner





SEM image of a silver contact pad, on the left in the overview, right a section of a pad corner

#### **SEM: FIB cut**

To avoid penetration of the measuring solution at the electrode edges, electrodes are designed to overlap the edge. In the example shown in Fig. 1, for example, this is  $3.8 \ \mu\text{m}$ . The extent of the overlap can be adjusted.





Left: SEM FIB cut of a gold electrode. Right: SEM FIB cut of a silver electrode